

#### **Better Crops with Arrays and Sequencing from Illumina**

#### AEIC, April 2013

Chandrasen Soans, Ph.D, MBA

Product Manager, Agrigenomics



#### **Focus Areas for Today's Presentation**

- Illumina-Our Vision
- The Toolset
- New Products and Services
- Trends and the Future



### Illumina at a Glance

- Founded in 1998
- ▶ IPO July 27, 2000
- Headquarters in San Diego, CA
- More than 850,000 sq. ft.
- Facilities in 7 countries
- Over 2,500 employees
  - ▶ >1200 R&D staff
  - ► >400 support personnel
- IP portfolio of >235 issued patents and 168 pending applications
- Added to the NASDAQ-100 listing in 2008
- Forbes "Fastest Growing Technology Company" 2007 & 2009





### We are a Global Organization

Expanded Manufacturing, R&D, Sales, Service & Support





#### **Vision Drives Change...Reality Exceeds Vision**



Scientists from the RAND Corporation have created this model to illustrate how a "home computer" could look like in the year 2004. However the needed technology will not be economically feasible for the average home. Also the scientists readily admit that the computer will require not yet invented technology to actually work, but 50 years from now scientific progress is expected to solve these problems. With teletype interface and the Fortran language, the computer will be easy to use.

5

# UNLOCKING THE POWER OF EVERY GENOME

#### Accelerating genetic gain through predictive breeding



To Targeted Validation and Beyond...

#### **Improving human health through agriculture**



# Having a Tangible Impact on Reduction Of Poverty And Hunger



- Improving rice productivity in South East Asia: IRRI
- Improving goat productivity in Africa: USDA & partners
- Developing Baobab tree genomics tools: Danforth Center, Monsanto & partners
- Improving Dairy productivity in Africa: University of New England (Australia)
- Improving productivity of pidgeonpea in India & Africa: ICRISAT
- Improving productivity of cassava in Africa: BecA

#### **We Have Seen an Unprecedented Publication Rate** 3739 peer-reviewed sequencing publications as of April 10, 2013





# **Enabling an Agrigenomics Revolution**

**Explosion Of Genomes Available At NCBI** 

http://www.ncbi.nlm.nih.gov/genome/browse/





# **Discover. Develop. Deploy.** *Access to multiple applications*



#### **Quick Portfolio Review**





### **Illumina Portfolio Overview**

#### **iNNOVATION** is in our DNA

#### From Genome-Wide Discovery to Targeted Validation and Screening



### HiSeq 2500





# **HiSeq 2500**

#### **Broadest Applications Flexibility**

Any application, any sample size. Operational efficiency.





#### MiSeq – Continuous Performance Improvements

Path towards 15Gb per run; enabling broader range of applications



"This information is intended to outline general product direction and it should not be relied on in making a purchasing decision. This material is for infogmation purposes only and may not be incorporated into any contract. This information is not a commitment, promise, or legal obligation to deliver this functionality. The development, release, and timing of any features or functionality described for our products remains at our sole discretion."

# MiSeq – R&D Demonstrated Scalability

#### Long read runs in excess of 20Gb



Output	22 Gb
Clusters	26.5M
Read length	450 (R1), 375 (R2)
Quality	71% ≥Q30
MR	1.7% (R1), 2.4% (R2)

Demonstrated by Illumina R&D



#### iScan<sup>TM</sup> System

#### POWERFUL & ECONOMICAL MICROARRAY PLATFORM



	Semi confocal laser scanning system			
<ul> <li>Resolution, 0.5 micron</li> </ul>				
Dual line - lasers at 532 & 658nm				
Fully automated operation				

Compact bench top system



# Add-on Content for Commercial Ag Genotyping Arrays

Product	Current BeadChip Format	Additional # Att BTs that can be supported on current BeadChip format	Maximum # of Beadtypes allowable on current substrate
BovineHD	8-sample	85,000	900,000
BovineSNP50v2	24-sample	31,000	90,000
CanineHD	12-sample	58,000	250,000
MaizeSNP50	24-sample	29,000	90,000
OvineSNP50	12-sample	183,00	250,000
PorcineSNP60v2	24-sample	32,000	90,000
BovineLD	24-sample	80,000	90,000
	Min # ABT: 1.0	00 Min # samples: 1.1	52



#### Steps in Building Genotyping Tools





# **Keeping Updated and/or Private Content**



# **Expanded Custom Genotyping Options**

Driving the Cycle of Discovery and Validation



Ability to add on supplemental content



# iSelect BeadChip Formats and Supported Marker Densities

BeadChip Format	24x1 2um	12x1 2um	4x1 2um
Number of Attempted Beadtypes (ABT)	3,072 – 90,000	90,001 - 250,000	250,001 – 1M
Add-on content range* (ABTs)	Up to limit of 90,000	Up to limit of 250,000	Up to limit of 1M
Instrumentation Needed		iScan, HiScan,	



# What is an Ag Consortium?



- "A consortium is an association of two or more individuals, companies, organizations or governments (or any combination of these entities) with the objective of participating in a common activity or pooling their resources for achieving a common goal."^
- A group of people or organizations with a common interest in advancing the collective understanding of the genetics of an agriculturally important species.
- The consortium members work together to enable actions that promote that interest

Consortia are helping advance genotyping in agriculture species around the world.

#### ^ Source: Wikipedia



# **Benefits Exist to Pooling Samples**

Economies of scale, also known as more for less!



When groups pool resources to create a genotyping tool, it lowers the barrier to entry by spreading the cost of synthesis of the beadpool over more samples. This was the case with the OvineSNP50 and BovineSNP50 commercial development projects. This was also the case with RosBREED Apple, Peach, Cherry and SolCAP Potato and Tomato genotyping beadchips.



25

Species	Identification	Classification	Provider	Consortium	# SNPs
Potato	Potato	Public	Illumina	SolCAP	8,303
Tomato	Tomato	Public	Illumina	SolCAP	7,720
Apple	Apple	Public	Illumina	RosBREED	8,788
Peach	Peach	Public	Illumina	RosBREED	8,144
Cherry	Cherry	Public	Illumina	RosBREED	5,696
Maize	MaizeSNP50	Commercial	Illumina	Commercial	56,110
Chicken	Chicken	Private: Public Sale	Illumina	Cobb Vantress- Hendrix-USDA	57,636
Cat	Feline	Private: Public Sale	Illumina	Morris Animal Foundation	62,897
Horse	Equine	Private: Public Sale	Illumina	Neogen (GeneSeek)	65,157
SheepHD	Ovine	Private: Public Sale	Illumina	AgResearch	680,000
Cattle	BovineHD	Commercial	Illumina	Various	777,962
Cattle	BovineSNP50v2	Commercial	Illumina	Various	54,609
Sheep	OvineSNP50	Commercial	Illumina	Various	52,241
Cattle	BovineLD	Commercial	Illumina	Various	6,909
Pig	PorcineSNP60	Commercial	Illumina	Various	62,163
Dog	CanineHD	Commercial	Illumina	Various	173,662
Wheat	Wheat 9k and 90K	Public	Illumina	Australia/US Wheat Consortium	8,632/81,587
Brassica	Brassica 50k	Public	Illumina	Intl Brassica SNP Consortium	52,157
Grape	Grape	Public	Illumina	ReSeq	18,071
Goat	Goat	Public	Illumina	Int'l Goat Genome	53,347

#### **New Products**





# IImoleculo

# Illumina Acquires Moleculo Inc. for Longer Reads

Pinit

January 8, 2013 by nextgenseek · 3 Comments



Illumina announced at the JP Morgan Healthcare Conference (JMP2013) that

it has acquired San

IIImoleculo

15

Share

Illumina Acquires Moleculo

Francisco based Mocleculo. Moleculo, the sequencing startup company (with not much information in the web), is a product of Stanford University's sequencing pioneer Stephen Quake.

#### Production of Long (1.5kb – 15.0kb), Accurate, DNA Sequencing Reads Using an Illumina

#### HiSeq2000 to Support De novo Genome Assembly

Date: Saturday, January 12, 2013 Time: 3:50 PM Room: Town and Country

Geoff Waldbieser, USDA-ARS, Stoneville, MS Michael Kertesz, Moleculo, Inc., San Francisco, CA Dmitry Pushkarev, Moleculo, Inc., San Francisco, CA Tim Blauwkamp, Moleculo, Inc., San Francisco, CA John Liu, Auburn University, Auburn, AL

Interspersed repeat sequences such as transposons and short tandem repeats shatter de novo genome sequence assemblies because short DNA sequences cannot span the repeat sequence. In order to produce a more contiguous genome assembly for the blue catfish, Ictalurus furcatus, we have used Moleculo's Long Reads product to generate extremely long and accurate reads through Illumina-based sequencing of libraries produced from long genomic fragments. To date, 201,508 long reads (933Mb total) have been produced from two libraries, ranged in length from 1.5kb to 15.8kb, and 83% of the total sequence was found in 129,554 long reads of at least 3.0kb. Pairwise alignments revealed 145,189 long reads contained only one or no mismatched bases along lengths of 400bp to 15,887bp. Preliminary assembly of only the long reads, using 99% sequence overlap identity, produced 46,098 contigs with an N50 length of 12.9kb and N80 length of 8.5kb. A further 42,141 long reads remained singlets (N80 = 4.6kb, N50 = 7.0kb). The long reads were aligned with the 1.6kb Tip1 and 1.0kb Tip2 transposons of channel catfish. Twenty one long reads contained an average 2.1kb (minimum 200bp) of sequence flanking the Tip1 orthologs. Similarly, 226 long reads contained an average 3.0kb of sequence flanking the Tip2 orthologs. The initial results demonstrate the utility of long and accurate DNA reads in bridging repetitive regions of the genome that cannot be otherwise resolved. Four additional Moleculo Long Read libraries are currently being sequenced to produce an additional ~2Gb of sequence in long reads.





# **Synthetic Long Reads**



#### Key components:

- Method for accurate amplification
- **Kit** for highly-parallel library prep
- Algorithm for long fragment reconstruction



#### NūPCR – Product Overview

**N**ūPCR (Gene Expression and Quantification)

Offerings: custom designed probe-based assays + Master Mix

- Illumina developed customizable assays and Master Mix
- DesignStudio custom primer and probe design



- qPCR DNA Binding Dye Assays (Gene Expression and Quantification)Offerings: custom designed primers
- Ready-to-use custom qPCR assays
- DesignStudio custom primer design

qPCR High Resolution Melt Assays (Genotyping)

Offerings: custom designed primers

- Ready-to-use custom HRM assays
- DesignStudio custom primer design





#### What is NūPCR?



- NūPCR is a probe-based qPCR technology that utilizes Illumina's unique NūZyme chemistry NūZyme is a multipart nucleic acid enzyme (DNAzyme).
- All oligos must assemble on the target sequence in order for fluorescence to be generated (high specificity). (5 oligos including primers)
- Universal substrate results in cost savings
- 4 dyes to choose fromallows multiplexing

### **Illumina Sample Prep Solutions**

Integrated workflows from sample to analyzed data



33

#### **Incorporation of dUTP in Second Strand Synthesis**



#### **Trends and the Future**





Sequence-Based Genotyping Next-Generation Genotyping Genotyping by Sequencing Genotyping using Next-Generation Sequencing



# **Sequence-Based Genotyping**

- Skim sequencing without reduced representation
- PCR and hybridization based methods of target enrichment
  - PCR based methods
  - Solid/Liquid phase hybridization methods
- Restriction enzyme methods of reducing representation
  - RADseq, GBS from Cornell and subsequent Cornell method



### **Genotyping by Sequencing** *Review of currently published methods*

Application Spotlight: Analyzing Genetic Variation

illumina

# Agrigenomics Genotyping Decisions Reach a Crossroads

For some applications, sequence-based genotyping provides a lower cost alternative to microarrays in performing genetic variation studies.

#### Introduction

Today's agrigenomics researchers have a wide variety of technologies at their disposal for collecting genetic information. Array-based methods of SNP screening have been the method of choice in analyzing and associating traits with regions of the genome for many plants and animals. As sequencing costs continue to drop, new approaches that leverage next-generation sequencing (NGS) technology are being developed to perform genotyping studies. We use the term next-generation sequencing-based genotyping (NGG) to encompass genotyping methods that leverage NGS technology. NGG includes targeted, reduced representation and hybridization-based approaches to discover and genotype SNPs, often simultaneously in many individuals or specimens. This application spotlight provides insight into different NGG methods, their benefits, and the role conventional array technology will play in the future.

#### Arrays Pave the Way in Agrigenomic Genotyping

In the late 1980s, researchers began identifying specific regions of DNA that influenced phenotypic traits in certain species. Efforts soon turned to the development of accurate and cost-effective genetic

populations<sup>4</sup>. By leveraging genetic screening, farmers and livestock breeders could gain immediate feedback, supporting better informed breeding decisions and accelerating their return on investment (ROI). Genotyping tools with a lower cost per sample could enable genetic screening to be performed routinely on large populations, with an attractive ROI offsetting the implementation cost of the technology.

#### Sequencing Advances Can Deliver More Cost-Effective Genotyping

The rapid evolution of sequencing technology has resulted in higher throughput and a lower cost per sample, often positioning NGG as a cost-effective and efficient agrigenomics tool for genotype screening, genetic mapping, purity testing, screening backcross lines, constructing haplotype maps, and performing association mapping and genomic selection<sup>5,8,7</sup>. The number of NGG methods continues to grow, with each offering the fundamental benefits sequencing provides, including reduced ascertainment bias, identification of variants other than SNPs (small insertions, deletions, and microsatelites), and an ability to perform comparative analysis across samples in the absence of a reference genome (Table 1).



# **Advantages of Sequencing-based Genotyping**

- Doesn't need a reference genome and can be used with an uncharacterized species
- Simultaneous identification and genotyping of SNPs
- Simple library prep needing low DNA input (100- 500 ng) and is generic
- High multiplexing lowers sequencing costs

#### **Points to Consider**

- Complex bioinformatics
- Missing data and therefore need to impute



#### Providing a Complete Solution with Arrays and Sequencing





Illumina is best positioned to provide you with the best tools, be it microarrays or sequencing



### **The Future**

- Accessibility to genomic tools will increase
- All economically impactful agricultural species, subspecies and their pathogens will be sequenced
- Field-based genome-wide diagnostics tools will be routine for pathogen and variety identification
- Genomic selection will surpass conventional methods as the dominant breeding paradigm





# **Questions?**

csoans@illumina.com



42