


Guidelines for Validation and Use of Immunoassays for Determination of Introduced Proteins in Biotechnology Enhanced Crops and Derived Food Ingredients

Dautlick
Grothaus
Hunst
Lipton
Magin
Mihaliak
Rubio
Stave



AEIC Analytical Environmental Immunochemical Consortium



About 27 member organizations, primarily:

- **agricultural biotechnology companies**
- **crop protection companies**
- **immunoassay kit manufacturers**

www.immunochem.org

AEIC Collaboration

- **Abraxis, LLC**
- **Dow AgroSciences, LLC**
- **Monsanto Company**
- **Pioneer Hi-Bred International**
- **Strategic Diagnostics Inc.**
- **Zeneca Ag Products**



Analytical Environmental Immunochemical Consortium




Mission:

To advance and support immunoanalytical technologies, the AEIC will:

- Provide a consensus voice for applications to agricultural, biotechnology and environmental industries
- Develop educational programs
- Furnish scientific expertise
- Establish performance standards for immunochemical methods

Immunoassay Evolution





Immunoassay

- **the specific binding of an antibody with its target analyte**
- **detection is achieved by the addition of a signal-generating component (reporter label)**
 - radioactivity (RIA)
 - enzymes (EIA, ELISA)
 - colloidal gold/latex
 - fluorescence
 - phosphorescence
 - chemiluminescence
 - bioluminescence




Sandwich ELISA



Capture Antibody



Analyte (protein of interest)



Labeled Antibody



Substrate




Product



Biotechnology Enhanced Crops

- **One or more genes are inserted into the plant genome by recombinant DNA techniques**
- **The new genes often code for the production of new proteins**
- **These introduced proteins (gene products) can be detected by immunoassay**



Ag Biotech Applications for Immunoassay

- **Ag Biotech Product Development**

- Screening/Transformant identification
- Line selection/Plant breeding
- Seed quality control (QC)



- **Regulatory Approval**

- Protein expression/Product characterization studies
- Assessments of food, feed and environmental characteristics
- Determining concentrations for toxicology studies




Ag Biotech Applications for Immunoassay

- **Approved Product**
 - Product support/stewardship
 - Intellectual property protection
 - Seed QC
- **Grain Handling/Distribution System**
 - Identification of transgenic products at grain elevator
 - Testing at other transfer points in supply chain
 - Identity preservation



Ag Biotech Applications for Immunoassay

- **Food Ingredient Testing**
 - Raw agricultural commodities
 - Processed food fractions



Method Validation and Use

- **Sampling**
 - Appropriate
 - Representative
- **Standardized Reference Materials**
 - Purified protein standards
 - Commodity to be tested containing known proportion of transgenic material



Method Performance Characteristics

- **Accuracy**
 - Recovery experiments
- **Extraction efficiency**
 - “Incurred” analyte vs. spiked
 - Exhaustive extraction
- **Precision**
 - Verify linearity of response
 - Measure variation among replicates



Method Performance Characteristics

- **Method Reproducibility**
 - Variation among replicates
 - within the test and among separate tests
 - among analysts
 - day to day
- **Sensitivity**
 - Limit of Detection (LOD)
 - response vs. background noise
 - Adequate sample size



Method Performance Characteristics

- **Specificity**
 - Interferences from reagents or labware
 - Matrix effects
 - Cross-reactivity
- **Ruggedness**
 - Independent validation
 - Round-robin or ring testing
 - Compare material from various sources




Considerations


- **Analyte stability**
- **Method sensitivity**
- **Difficult matrices**
- **Defined testing objective**
 - Qualitative or yes/no
 - Semi-quantitative or threshold
 - Quantitative




Microplate ELISA




STANDARD CURVE FOR SANDWICH ELISA



Immunoassay Strip Test Lateral Flow Device



Lateral Flow Strip Immunoassay



Conclusion

- **Immunoassays are extremely useful for Ag Biotech**
 - both microplate and lateral flow devices
 - for quantitative and threshold applications
 - for product development and support
- **Like any analytical method**
 - Immunoassays must be validated
 - Performance characteristics demonstrated and described



Acknowledgements

- AEIC
- Leah Porter
- Glen Rogan
- Dean Layton
- Jim Rittenburg
- Mike Morgan and the Journal of Food and Agricultural Immunology

